ENGLISH

INSTRUCTIONS FOR USE

1. DESCRIPTION

BD PosiFlush™ XS Syringe is a ready to use sterile medical device (according to 93/42/EEC medical device directive). It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. The contents of our unopened or undamaged blister packages is guaranteed to be sterile.

2. INTENDED USE

- BD PosiFlush™ XS Syringe is intended **FOR FLUSHING ONLY** in-situ vascular access devices.
- BD PosiFlush™ XS Syringe is not intended for dry product reconstitution, for medication dilution, or where intravenous therapy with sodium chloride is indicated.
- Using aseptic technique, **BD PosiFlush™ XS Syringe can be used on a sterile field.**

3. CONTRAINDICATIONS

None known.

4. CAUTIONS

Do not re-sterilise before use.

Do not use if unit package or content is damaged. Verify the expiration date on the product package or label. Do not use if product has expired.

Do not use if syringe tip cap or stopper is damaged in any way that suggests or indicates syringe leakage.

Do not use if solution contains a precipitate, or has any type of suspended particulate matter.

For single use only.

Discard any partially used product.

Store at controlled temperature (15-25°C). Excursions permitted to 30°C.

Do not leave at freezing temperature.

Federal (USA) Law restricts this device to sale by or on the order of a physician.

Re-use may lead to infection or other illness/injury.

Check with drug manufacturer instructions for use to ensure compatibility with 0.9% sodium chloride solution prior to use. If 0.9% sodium chloride solution is not compatible, follow the drug manufacturer instructions for flushing practices, or first flush the vascular access device (VAD) with a compatible solution such as 5% dextrose in water to remove traces of the medication followed by flushing with 0.9% sodium chloride solution as a locking solution.

Clinicians should consider the patient's specific medical conditions, treatment needs, age, and weight that may

require restricted sodium or fluid intake when deciding to flush with 0.9% sodium chloride Injection.

5. DIRECTIONS FOR USE

Single use, single patient device. To ensure safe medication preparation and administration, clinicians should practice the "7 rights" of medication administration: right patient, right drug, right dose, right time, right route, right reason and right documentation.

Use aseptic technique throughout the procedure.

1. *Open blister pack and remove syringe.*

Place syringe on a sterile field if applicable.

2. *Check that syringe tip cap is in place. Inspect clarity of solution. (Fig. 1a)*

3. *Depress plunger with tip cap on to release the stopper seal. (Fig. 1b)*

4. *Unscrew tip cap from the syringe ensuring that there is no touch contamination of the syringe luer connection. (Fig. 2)*

5. *Push syringe plunger to expel the air. (Fig. 3)*

6. *Connect BD Luer-Lok™ syringe to vascular access device, taking care that there is no touch contamination of the connection. (Fig. 4)*

7. *Push syringe plunger to flush the required volume of saline following institution's policy. (Fig. 5)*

In case of extreme plunger resistance, it is recommended that excessive force is not exerted. Some patients may experience a transitory taste or odour during flushing. This minor effect ceases shortly after the procedure is completed.

8. *After use, dispose of in accordance with recognised procedure in your institution. (Fig. 6)*

DO NOT RE-USE.

6. COMPOSITION PER UNIT

Polypropylene syringe with BD Luer-Lok™ tip. Elastomeric stopper free of natural rubber latex. Saline solution: sodium chloride 9g/l (NaCl 0.9%), distilled sterile water to volumes, preservative free and non-pyrogenic.

Single blister pack consists of paper and Polypropylene - Polyamide film.

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BD, the BD Logo and PosiFlush are trademarks of DGP4/904
© 2020 BD. All rights reserved. Rev. 2020-01
EC REP Becton, Dickinson & Co. Ltd. Donore Road, Drogheda, Co.Louth, Ireland