Catheter Stabilization Device

Read carefully before use

Safety and Efficacy Considerations:

- Single use only.
- Do not re-sterilize.

Rx Only, Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

STERILE Sterile by ethylene oxide. Do not use if package is damaged. Fluid path (if present) is nonpyrogenic.

REF Product code, **LOT** lot number and expiration date information – see packaging.

Do not alter the StatLock® stabilization device or components. Procedure must be performed by trained personnel with knowledge of anatomical landmarks, safe technique and potential complications.

LATEX-FREE

Indications for Use:

The StatLock® stabilization device is for compatible medical tubes and catheters. Multiple pad and retainer designs are available.

Choose the StatLock® stabilization device indicated for the tube or catheter to be secured.

Contraindications:

Known tape or adhesive allergies.

Warnings and Precautions:

1. Do not use the StatLock® stabilization device where loss of adherence could occur, such as with a confused patient, diaphoretic or nonadherent skin, or when the access device is not monitored daily.
2. Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the StatLock® stabilization device.
3. SUTURE THE StatLock® STABILIZATION DEVICE PAD TO SKIN IF DESIRED OR DEEMED NECESSARY.
4. Avoid contact with alcohol or acetone, both can weaken bonding of components and pad adherence.
5. Minimize catheter/tube manipulation during application and removal of the StatLock® stabilization device.
6. A StatLock® stabilization device Luer-Lock connector must be used to secure central venous and arterial catheters.
7. Always apply adhesive strip to central venous and arterial catheters at or near insertion site when using StatLock® stabilization device.
8. Remove oil and moisturizer from targeted skin area.
9. StatLock® stabilization device pad adherence and catheter/tube position should be routinely inspected.
10. Orient StatLock® stabilization device so arrow points toward catheter tip.
11. StatLock® stabilization device should be replaced at least every 7 days.
12. This is a single use device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
13. Do not re-sterilize. The sterility of the single use device is not guaranteed following re-sterilization because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Re-sterilization may compromise the structural integrity, essential material and/or design characteristics and may lead to an unpredictable loss of functionality and/or device failure.

Use sterile technique for StatLock® stabilization device application and removal. See illustrations for specific details.

Extension Set: If present, Prime/Flush/Connect per accompanying instructions.

StatLock® stabilization Device Application:

1. **PREP:** Prepare targeted area using alcohol, to remove oil and moisturizer, and skin prep for skin protection and enhanced adhesion. ALLOW TO DRY COMPLETELY AS CONTACT WITH ALCOHOL CAN WEAKEN BONDING OF COMPONENTS AND PAD ADHERENCE. ALWAYS CONNECT THE StatLock® STABILIZATION DEVICE TO CATHETER BEFORE ADHERING PAD TO SKIN.
2. **PRESS:** Before applying the StatLock® stabilization device pad to skin, secure tube or catheter to StatLock® retainer.
3. **PEEL:** Peel off the StatLock® stabilization device paper backing, and
4. **PLACE:** Place the StatLock® stabilization device on targeted skin one side at a time.

Note:

- The StatLock® stabilization device pad can be adhered after X-ray confirmation of proper catheter/tube position.
- Apply adhesive strip to central venous and arterial catheters for added security during StatLock® stabilization device dressing changes.

The StatLock® stabilization device general removal:


Carefully remove overlying dressing using stretch technique and sutures (if applicable).

1. **DISENGAGE** catheter/tube from StatLock® stabilization device retainer.
2. **DISSOLVE** StatLock® stabilization device adhesive with alcohol swab while gently lifting StatLock® stabilization device pad.
3. **DOCUMENT** StatLock® stabilization device dressing change in the patient chart.
4. **DISPOSE** of all equipment in appropriate containers.

Extension Set

Indications for Use:

For use as an extension to the IV fluid delivery system. StatLock® stabilization device extension set male Luer adapter is specifically designed to secure the extension set when used as part of the StatLock® stabilization device system.

Directions for Use:

Use aseptic technique.

1. Tighten all connections.
2. Remove Protector Caps (A) as required.
3. Connect Female Port (B) to fluid delivery line or device and prime/flush.
4. Connect Male Port (C) to catheter hub.
5. Open Slide Clamp (D) for unimpeded flow.
6. Swab septum of injection site and SmartSite™ device (if present) with preferred antiseptic.

General Precautions:

1. Avoid air trapped in set. Puncturing tube can cause air embolism or leakage.
2. Avoid contact with alcohol, as contact with alcohol can weaken bonding of components and pad adherence.

Options / Precautions:

Needless Valve (E) and Y-Site w/ Needless Valve (F)

1. Connect fluid line or device to center of valve using Luer Lock connections only.
2. Do not use with a blunt tip cannula or conventional needle. If a needle must be used, attach a PRN adapter and use as a standard injection site.

T-Connector (G)

Stopcock (H)

1. Position handle over port you wish to close.
2. Open port for fluid delivery.
3. Close and cap available port when not in use.

Standard Y-Site (I)

1. Use with a conventional needle only. Do not use with a blunt tip cannula.

NOTES

1. Tubing is non-DEHP PVC
2. Dispose used set in appropriate container.
3. Replace extension set and ports per CDC guidelines.

Micro	Fill Volume (mL)	
	No Valve	W/Valve
12" (304.8mm) ext set	0.40	NA
4" (101.6mm) ext set	0.18	NA
6" (152.4mm) ext set with SmartSite®	0.25	0.34
6" (152.4mm) ext set with T-Connector and SmartSite®	0.18	0.27
8" (203.2mm) ext set with 3-way stopcock	0.36	0.55
8" (203.2mm) ext set with bonded 3-way stopcock	NA	0.55
9" (228.6mm) bifurcated ext set with SmartSite®	0.57	0.75

Macro

6" (152.4mm) ext set with Smart Y® and SmartSite®	0.90	1.06
6" (152.4mm) ext set with SmartSite®	0.90	0.99

Bard and StatLock are trademarks and/or registered trademarks of C. R. Bard, Inc.

All other trademarks are the property of their respective owners. © 2014 C. R. Bard, Inc. All Rights Reserved.

Manufacturer:
Bard Access Systems, Inc.
605 North 5600 West
Salt Lake City, UT 84116 USA
800-443-5505, 801-522-5000
www.bardaccess.com
Assembled in Mexico

Authorized Representative in the European Community
Bard Limited
Forest House, Brighton Road
Crawley West Sussex
RH11 9BP UK

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Bard et StatLock sont des marques et (ou) des marques déposées de C. R. Bard, Inc.
© 2014 C. R. Bard, Inc. Tous droits réservés.