# Umbilical cord clamp

Clamp para cordão umbilical

ZILLION<sup>®</sup>  
MEDICAL DISPOSABLES




Bastos Viegas, s.a.

IFU55199910. Rev. 2 (03-2025)

Avenida da Fábrica, 298, 4560-164 Guifufe, Penafiel - Portugal

Tel.: 255 729 500 - geral@bastosvegas.com - www.bastosvegas.com

<https://www.bastosvegas.com/downloads/fiu/IFU55199910.pdf>

## EN UMBILICAL CORD CLAMP

INTENDED PURPOSE, INDICATIONS AND CLINICAL BENEFIT

Umbilical cord clamp is intended to clamp the umbilical cord after delivery to prevent bleeding from the cord stump after the cord had been cut.

### INSTRUCTIONS FOR STERILIZATION

The device intended to be sterilized must be packaged within a sterile barrier system (SBS) permeable and compatible with the sterilization parameters that will be described below. The system shall comply with EN ISO 11607-1/2 and CEN ISO/TS 16757:2021. It is the responsibility of the processor to select an appropriate SBS and carry out the sterilization process.

The following EO sterilization process parameters have been validated according to EN ISO 1135, for:

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CYCLE PARAMETERS SPECIFICATION DURATION

PRECONDITIONING

Minimum temperature during the preconditioning 21,4 °C

Temperature 33 - 46 °C

Humidity 50 ± 15 %RH

STERILIZATION

50 mbar

Humidification/Dwell 60 - 95 %RH

Gas injection 50 mbar of N<sub>2</sub>

Exposure pressure 930 mbar

Exposure EO concentration 580 - 750 mg/L

EO removal temperature 48 - 55 °C

Flushing 100 mbar

Temperature 70mbar x 3cycles

AERATION

70 - 20 m/s

Speed of air extraction

Two cycles with the above parameters have been validated. The second sterilization cycle is only allowed for unused devices. After undergoing sterilization processing, the device's shelf life will be limited to 5 years. This 5-year durability corresponds to the time the device is stored in the sterile environment and the validation of the non-sterile component. EO sterilization processing shall be performed by trained personnel and by utilizing adequate and validated EO equipment. The processor shall comply with EN ISO 1135 and any national legislation where the processing will occur. If transport is applicable, it must be validated acc. to international applicable standards (e.g. ASTM4169, ISTA).

### INSTRUCTIONS FOR USE

After sterilizing:

- 1-Check the thickness of the umbilical cord. In case it is a thin cord, it may be necessary to fold it before placing the clamp to obtain a bigger resistance to the clamp closing.
- 2-Place the cord within the clamp, in the clamping zone (figure 1a).
- 3-Check that there is not any object blocking the locking site (e.g. swabs, surgical instruments, etc.).
- 4-Close the clamp until it locks and clicks (figure 2). Make sure that the clamp is well closed and that it is impossible to open it with the hands. Cut the umbilical cord.
- 5-Confirm that there is no bleeding at the end of the cord. Check regularly for signs of bleeding due to a misplaced clamp or due to a damaged cord.

### WARNINGS AND PRECAUTIONS

• Intended to be used exclusively by paediatrics/ neonatology healthcare professionals.

• Visually inspect the device. Do not use the device if you detect any sign of damage or defect.

• Always check the integrity of the locking system before use. Do not use if incomplete or blocked. Incorrect placement of the clamp on the umbilical cord may lead to an inadequate seal and consequent blood loss that might prove fatal for some babies.

• Do not use the device if it is used for a medical purpose other than its intended use.

• When used on healthcare facilities, and if contaminated:

• Dispose according to applicable legislation and/or infection prevention standards for human waste, to minimize potential risk of contamination, infection or pollution.

• Otherwise: Dispose according to applicable legislation and/or infection prevention standards for urban waste.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents can occur with the use of devices that are not sterilizing or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your country.

LABELLING SYMBOLS

1)UKCA marking 2)United Kingdom authorized representative

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