Steriflush®

Prefilled Syringes

Steriflush Medical Products
Zona Industrial 1 - Lote 11 a 14
4560-164 Gualtufe, Penafiel
PORTUGAL www.steriflush.com

CE 0344

STI

EN

ENGLISH - For flushing only
Prefilled syringes 0.9% Sodium Chloride - Saline

It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. STERIFLUSH® prefilled syringe 0.9% Sodium Chloride is intended for FLUSHING ONLY in-situ vascular access devices.

Cautions:

- STERIFLUSH® prefilled syringe 0.9% Sodium Chloride is not intended for reconstitution of drugs or patient rehydration, for medication dilution, or where intravenous therapy with sodium chloride is indicated.
- Do not use if solution is discoloured.
- Do not use if the solution contains a precipitate, or has any type of suspended particulate matter.
- Do not use without carefully inspecting and discard if does not meet visual specifications.
- Do not use if syringe tip cap or stopper is damaged in any way that suggests or indicates syringe leakage.
- Some patients may experience a transitory taste or odour during flushing. This minor effect ceases shortly after the procedure is completed.
- Once opened, the Medical Device must be used as soon as possible.

Instructions for use:

- STERIFLUSH® prefilled syringe 0.9% Sodium Chloride must be used in line with your hospital's intravenous drug administration and vascular access policy and also the vascular access device manufacturer's guidelines.
- Using aseptic technique, remove the syringe from the package. Discard if packaging is not intact.
 - Check that syringe tip cap is in place. Inspect clarity of solution (Fig. 1).
 - Depress plunger with tip cap on to release the stopper seal (Fig. 2).
 - Unscrew tip cap from the syringe ensuring that there is no touch contamination of the syringe luer connection (Fig. 3).
 - Push syringe plunger to expel the air (Fig. 4).
 - Connect carefully STERIFLUSH® syringe to vascular access device, taking care that there is no touch contamination of the connection (Fig. 5).
 - Push syringe plunger to flush the required volume of saline following institution's policy (Fig. 6). In case of extreme plunger resistance, it is recommended that excessive force is not exerted, as this may damage the syringe.
 - After use, dispose of in accordance with recognised procedure in your institution (Fig. 7).

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