Lipifil Air 24 h, PVC-free

REF 2909301

de Gebrauchsinformation

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en Instruction for use

Please read this information carefully before use

Product name

Lipifil Air 24 h infusion filter with integrated PET membrane, pore size 1.2 µm (indicated by yellow arrows), integrated PTFE venting membrane, pore size 0.03 µm and tube connection with Luer-Lock connectors.

Intended use

The filtration of lipid emulsions and compound solutions to retain micro-organisms including fungi, spores, and particles (>1.2 µm) and to separate air bubbles. Commercially available infusion sets and connecting tubes with Luer-Lock hubs are used.

Duration of use

Depends on intended use up to a maximum duration of 24 hours. Blocked filters must not be rinsed, as the membrane could be damaged and existing contamination flushed through. The blocked filter must be removed and replaced. Filter blockage is an indication of contamination or incompatibility of medications. Injection of medications can result in excessive pressure on the filter surface and destruction thereof. In this case, please use special filter systems with an injection port.

Restrictions on use

The infusion filter must only be used as intended for infusion therapy. Blood, cellular components, coagulation factors, and suspensions must not be administered via this filter.

Warning notes

Do not use product after the expiry date printed on the package.

This product does not replace necessary hygienic and aseptic measures.

Take precautions to prevent infection.

Integration into an infusion set

1. Remove filter from sterile packaging.
2. Remove proximal protective cap.
3. Connect vented administration set to filter inlet.
4. Remove distal protective cap.
5. Connect infusion line to filter outlet.
6. To vent filter with infusion line:
 - Hold up the free connection of the distal end of the filter
 - Fill system with connected infusion solution

- until this emerges at end of infusion line, filling filter from bottom to top for proper venting.
7. Connect distal end of the filter to the patient’s IV access.
 8. Check integrity of system (visual check).
 9. Dispose of protective caps, packaging and other materials no longer needed.
 10. Start infusion.
- Materials used**
Flat infusion filter housing and packaging
PVC-free flat infusion filter housing; transparent MABS
PET filter membrane.
Protective caps and sealing cones: PE and ABS
Tube connections: PUR
Dispose according to hospital protocol.

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Fresenius Kabi AG