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Tel: 0212 622 04 00

www.setmedikal.com.tr

GB CATHETER TIP HYPODERMIC SYRINGE INSTRUCTION FOR USE

Catheter Tip Hypodermic Syringe
Catheter 3 part Syringe 50ml/ 60 ml

Product Type	Catheter Tip Hypodermic Syringe
Catheter – Centric 3 parts	

Intended Use

Hypodermic Syringes, made of empty, sterile, and disposable plastic materials, are used by end users for the aspiration and injection of liquids after filling. The syringes are primarily for human use. Sterile syringes are designed to be used immediately after filling and are not intended to retain medication for long period of time.

Sterile, single-use hypodermic syringes with catheter tips are devices used for the administration of medical products and fluids.

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PP – Polypropylene
Gasket	Isoprene rubber, Latex free.
Cap	PP – Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Instruction for Use

- Open the package by peeling it from the place marked by the arrow.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging includes the sterile syringe.
- Make sure that the grading lines on the syringe cylinder are legible.
- Determine the amount of liquid to be used using the graduated scale on the syringe.
- Make sure the syringe is free of air bubbles. If there is any air bubble in the syringe, be sure to remove it.
- Before injection, check if fluid is flowing in the syringe.
- The electronic version of PR.A112 Instruction for use is available at www.setmedikal.com.tr.
- The capacity tolerance of the graduated scale on the syringe is stated in the table below in accordance with EN ISO 7886-1 Standard.

Nominal capacity of syringe V*, (ml)	Graduated Capacity Tolerance (%)
5≤V≤60	±4

V*: Volume

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Device Description Intended for Use With The Device

It is used with medical devices that will be flushed, such as Foley catheters and gastrostomy tubes.

Indications

Catheter tip syringes are devices used to administer medical products and liquids. Syringes are designed to be used for intravenous, intramuscular, subcutaneous and intradermal injection, intravenous blood collection and feeding medical systems. The product shelf life is 5 years.

Contraindications

Catheter tip syringes should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or resterilization is not possible. It is not intended for use other than its stated intended use. Dispose of as medical waste after use.

Cautions

- Do not use the Product if the package is opened and damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- When the sterile barrier system (unit packaging) is opened, the needle is inserted without waiting in accordance with aseptic rules.
- If there is a special warning for the drug administered by syringe, it should be considered during application.
- In applications, the capacity tolerance table of the graduated scale given in the instructions for use should be considered.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination, and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Responsible Person in the United Kingdom		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		Marking for meeting the requirements of UK Medical Device Regulation The number displayed on the bottom is the Notified Body Number
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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