Umbilical cord clamp

Clamp para cordão umbilical


- FR** Clamp de Bar umbilical
- DE** Nabelschurklemme
- NL** Navelstrengklemme
- IT** Morsetto per cordone ombelicale
- ES** Pinza para cordón umbilical
- SV** Navelsträngsklämma
- GR** Κλάμμα za πλῆνα βρῆς
- HR** Stezaljka za pupčanu vrgcu
- SI** Stovzla na pupečnik
- DA** Navelstrængsklemmer
- EN** Nabanöör klamber
- EL** Λαπαρόνηο ρυτίστη
- EL** Λαβίδα ομφάλτου λώρου
- EL** Κολόκλιζιόνσπρεζ
- LV** Nabassaites spalles
- LT** Virkėsties spaustukas
- PL** Zakęś do pęrowiny
- RO** Cleamă cordon ombilical
- SK** Svorka na pupečník šnúru
- SL** Spjonka za porokovino

https://www.youtube.com/watch?v=SXUVjPxaUBQa

Bastos Viegas, s.a.

Avenida da Fábrica, 298, 4560-164 Guilhufe, Penafiel - Portugal
Tel.: 255 729 500 - geral@bastosviegas.com - www.bastosviegas.com

IMPORTANT

CAUTION: incorrect placement of the clamp on the umbilical cord may lead to an inadequate seal and consequent blood loss that might prove fatal for some babies. Read the instructions before using the device.

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EN UMBILICAL CORD CLAMP

INTENDED PURPOSE, INDICATIONS AND CLINICAL BENEFIT

Umbilical cord clamp is intended to clamp the umbilical cord after delivery to prevent bleeding from the cord stump after the cord had been cut.

INSTRUCTIONS FOR STERILIZATION

The device intended to be sterilized must be packaged within a sterile barrier system (SBS), permeable and compatible with ethylene oxide (EO) sterilization and with the sterilization parameters that will be described below. The system shall comply with EN ISO 11607-1/2 and CEN ISO 15765:2021. It is the responsibility of the manufacturer to select an appropriate method of packaging under controlled environmental conditions. The following sterilization process parameters have been validated according to EN ISO 11335, for:
-Density sterilization load of 180 kg/m³ (in euro pallets);
-Use of a headcap bag and a paper/film peel-pack and card and carton box;
-Header bag SBS (acc. to EN ISO 11607-1/2 and EN 868-5); and
-Carton and carton box protective packaging,
-Removal of the SBS (acc. to EN ISO 11607-1/2 and EN 868-5), and
-Flushing.

CYCLE PARAMETERS	SPECIFICATION	DURATION
Minimum temperature at preconditioning	21.4 °C	-
Temperature	33 - 46 °C	22-72h
Humidity	50 ± 15 %RH	-
STERILIZATION		
Vacuum/Air removal	60 - 95 %RH	60min.
Humidification/Dwell	350 mbar of N, 480 mbar of EO	52min.
Gas Injection	930 mbar of N	-
Exposure pressure	580 - 750 mpa/L	420-1335min.
Exposure EO concentration	48 - 55 °C	-
Exposure temperature	100 mbar	17min.
EO removal	70mbar x 3cycles	-
Flushing	ALÉRIATION	-
Temperature	28 - 42 °C	Min.1th
Speed of air extraction	10 - 20 m/s	-

Two cycles with the above parameters will be performed. The second sterilization cycle is only allowed for the sterilization of the device. After sterilization processing, the device's shelf life will be limited to 5-years. This 5-year duration corresponds to the shorter shelf life between the sterilized component and the validity of the non-sterile component. EO sterilization processing should be performed by trained personnel and by utilizing adequate and validated EO equipment's. The processor shall comply with EN ISO 11335 and any national or local legislation where the processing will occur. If transport is applicable, it must be validated acc. to international applicable standards (e.g., ASTM4169; ISTA).

INSTRUCTIONS FOR USE

- After sterilizing:
 - 1-Check the thickness of the umbilical cord. In case it is a thin cord, it may be necessary to hold it before placing the clamp to obtain a bigger resistance to the clamp closing.
- Place the cord within the clamp, in the clamping zone (figure 1a).
- Check that the device is not in any object blocking the locking site (e.g. swabs, surgical instruments, etc.).
- Close the clamp until it locks and clicks (figure 2). Make sure that the clamp is well closed and it is possible to open with the hands. Cut the umbilical cord.
- Confirm that there is no bleeding at the end of the cord. Check regularly for signs of bleeding due to a misplaced clamp or due to a damaged cord.

WARNINGS AND PRECAUTIONS

- Intended to be used exclusively by paediatric/s neonatology healthcare professionals.
- Intended to be used sterile.
- Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Always check the integrity of the locking system before use. Do not use if incomplete or if the locking system is damaged. This may lead to an inadequate seal and consequent blood loss that might prove fatal for some babies.
- Do not re-use. Intended to be used one single time to prevent risks of contamination, including possible infections, and/or functional damage.
- When used on healthcare facilities, and if contaminated:**
 - Dispose according to applicable legislation and/or prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.
- Otherwise:**
 - Dispose according to applicable legislation and/or infection prevention standards for urban waste.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise essential safety or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

LABELLING SYMBOLS

1)UKCA marking 2)United Kingdom authorized representative

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