DASH 6™ Syringes (BNS)

INSTRUCTIONS FOR USE

	Page
EN DASH 6™ ISO 80369-6 NRFit® Compliant Syringe	4
FR Seringue DASH.6™ conforme à la norme ISO 80369-6 NRFit®	5
IT Siringa DASH 6™ conforme a ISO 80369-6 NRFit®	6
ES Jeringa NRFit® DASH 6™ que cumple la norma ISO 80369-6	7
DE DASH 6™ ISO 80369-6 NRFit®-konforme Spritze	8
DA DASH 6™ ISO 80369-6 NRFit®-kompatibel sprøjte	9
NO DASH 6™ ISO 80369-6 NRFit®-kompatibel sprøyte	10
SV DASH 6™ ISO 80369-6 NRFit®-kompatibel spruta	11
NL DASH 6™-NRFit®-sput die voldoet aan ISO 80369-6	12
PL Strzykawka DASH 6™ zgodna z NRFit® ISO 80369-6	13
PT Seringa DASH 6™ NRFit® compatível com a norma ISO 80369-6	14
EL ΣΥΡΙΓΓΕΣ DASH 6™ ΣΥΜΒΑΤΕΣ ΜΕ ΤΟ ISO 80369-6 NRFit®	15
TR DASH 6™ ISO 80369-6 NRFit® Uyumlu Şiringa	16
CS DASH 6™ – Stříkačka podle normy ISO 80369-6 NRFit®	17
ET DASH 6™ ISO 80369-6 NRFit®-ühilduvusega süstal	18
HU DASH 6™ ISO 80369-6 NRFit® kompatibilis fecskendő	19
LV ISO 80369-6 standartam atbilstoša NRFit® šļirce DASH 6™	20
MT Siringa DASH 6™ Konformi mal-ISO 80369-6 NRFit®	21
RO Seringă DASH 6™ compatibilă ISO 80369-6 NRFit®	22
SL BRIZGALKA DASH 6™, skladna s standardom ISO 80369-6 NRFit®	23
SK Injekčná striekačka DASH 6™ ISO 80369-6 NRFit®	24
BG Спринцовка DASH 6™, отговаряща на стандарта ISO 80369-6 NRFit®	25
LT DASH 6™ standartą ISO 80369-6 NRFit® atitinkantis švirkštas	26

DASH 6™ ISO 80369-6 NRFit®**Compliant Syringe****Instructions For Use**

Please read these instructions before use.

INTENDED USE

The DASH 6™ NRFit® Lock Syringe and DASH 6™ NRFit® Slip Syringes are intended to be used with ISO 80369-6 NRFit® neuraxial compliant devices for the administration of neuraxial medication or anaesthetic.

INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial device and is used to administrate neuraxial medication and anaesthetic to a patient.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for intravenous or enteral use.

WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- Do not use if there are no dose markings on the syringe.
- Do not re-use or reprocess the device as this may lead to leakage or infection.
- Do not use if the tip is visibly blocked or the device is damaged.
- The device is single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device.

- Verify that the neuraxial device being used is compatible with DASH 6™ ISO 80369-6 NRFit® to reduce the risk of misconnections.
- Fully depress the plunger before filling with neuraxial medication/ anaesthetic. After filling, ensure that the syringe moat is free of excess neuraxial medication/anaesthetic and that there are no air bubbles in the syringe. If there is excess neuraxial medication/ anaesthetic or there are any air bubbles, remove them by flicking the side of the syringe.
- Attach the ISO 80369-6 NRFit® compliant syringe into the compatible medication port.
- Once this is secure, gently push the syringe plunger to transfer the medication/anaesthetic into the neuraxial site.
- After completely expelling the medication/anaesthetic from the syringe, detach the syringe and dispose of immediately as per hospital protocol.
- It is recommended to use a compatible filter in conjunction with the neuraxial device, or a drawing up device with an integrated filter, to prevent entrainment of any foreign particulate matter when drawing up fluid into the device.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

GRADUATION ACCURACY

Graduation accuracy detailed below complies with the requirements of ISO 7886-1.

Syringe volume (Nominal Capacity)	Capacity tolerance (Greater than or half-filled syringe)	Capacity tolerance (Less than half-filled syringe)
1ml	± 5% of expelled volume	± 1.5% of nominal syringe capacity x ± 2% of expelled volume
3ml	± 5% of expelled volume	± 1.5% of nominal syringe capacity x ± 2% of expelled volume
5ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
10ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
20ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
60ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume

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 **GBUK Group Ltd.**
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukgroup.com
info@gbukgroup.com

EU REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland

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Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukgroup.com.

Not made with natural rubber latex.

DASH 6™ is a trade mark of GBUK Group Ltd., registered in the UK and other countries and regions.

NRFit® is a registered trade mark of GEDSA and used with their permission.

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**GBUK GmbH**