FR	Instructions d'utilisation des robinets ou rampes avec et sans prolongateurs en PE/PVC et PVC
EN	Instructions for use of stopcocks or manifolds with and without PE/PVC and PVC extension tubing
ES	Instrucciones de uso de las llaves o rampas con y sin prolongadores de PE/PVC y PVC
IT	Istruzioni per l'uso dei rubinetti o delle rampe con o senza prolunghe in PE/PVC e PVC
PT	Instruções para uso de torneiras ou rampas com e sem prolongadores de PE/PVC e PVC
DA	Brugsvejledning for haner eller ramper med og uden forlængere i PE/PVC og PVC
NL	Gebruiksaanwijzingen voor kranen of kranenbalken met en zonder verlengslangen van PE/PVC en PVC
SL	Navodila za uporabo ventilov in razdelilnikov z ali brez podaljškov iz PE/PVC in PVC
SV	Instruktioner för användning av kranar eller ramper med och utan förlängningsslangar PE/PVC och PVC
HR	Upute za uporabu skretnica ili rampi s ili bez nastavaka izrađenih od PE/PVC-a i PVC-a
FI	Jatkoliittimillä varustettujen ja varustamattomien venttiilien tai kokoojaputkien käyttöohjeet, materiaali PE/PVC ja PVC
TR	PE/PVC ve PVC uzatma hatları bulunan ve bulunmayan musluklar veya çıkış boruları için kullanım talimatları
DE	Bedienungsanleitung für Hähne oder Gestelle mit und ohne Verlängerungen aus PE/PVC und PVC
NO	Instrukser for bruk av ventiler og broer med eller uten forlengerslanger i PE/PVC og PVC
PL	Sposób użycia zaworów lub listew z przedłużaczami i bez nich z PE/PVC i PVC
HU	Szelepek és elosztók használati utasítása PE/PVC és PVC hosszabbítókkal vagy azok nélkül
EL	Οδηγίες χρήσης στροφιγγών ή ραμπών με ή χωρίς επεκτάσεις από PE/PVC και PVC
AR	تعليمات استخدام صناعير أو أذرع الرش بدون أو بوصلات إطالة مصنوعة من البولي إيثيلين/البولي فينيل كلوريد ومن البولي فينيل كلوريد (الرقم المرجعي:
CS	Pokyny k použití kohoutků či ramp s prodlužovacími hadičkami či bez nich z PE/PVC a PVC
SK	Návod na použitie ventilov alebo rámp s nadstavcami PE/PVC a PVC

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1. Indications/Use claims:

3-way stopcocks and manifolds are devices intended to be placed on an infusion line. They allow administration of multiple drug substances simultaneously or in alternation.

These devices are intended for use by healthcare professionals only.

These devices are intended for use on pediatric or adult patients.

2. Composition:

Stopcocks or manifolds with and without extension tubing may include non-return valves or bidirectional valves.

Extension tubing with stopcocks or stopcock manifolds can have different diameters, lengths or materials (please refer to technical data sheets).

These devices do not contain DEHP (DOP), latex, BBP, DBP, DIBP or DHP. Any presence of these compounds is a trace amount and is compliant with REACH regulations.

3. Clinical benefit/Performance/Action mechanism:


Stopcocks and manifolds for intravenous infusion are sterile, single-use medical devices MDs.

Action mechanism:

1: the stopcock is closed on one side (left on the diagram) and open on the other two sides, allowing the drug substance to flow in both these directions.

2: the stopcock is closed on one side (bottom on the diagram) and open on the other two sides, allowing the drug solution to flow in a horizontal direction.

3: the stopcock is open on all three ports.

A : non-watertight stopper, must be removed and discarded before connecting the patient to the infusion line

4. Prerequisites before use and instruction for use:

Comply with applicable hygiene rules.

Check the integrity and cleanliness of the packaging and the device before use, in case of defect do not use the device.

Gently open the blister, taking into account the direction of peeling. Check the expiration date before use, if expired do not use.

Check all connections and stopper caps fully tightened before using the device. Purge the entire device before connection to the patient. If non-watertight (coloured) stoppers are used, they must be removed and discarded before connecting the patient to the infusion line. Purge stopcock after stopcock, all ports, beginning with the distal end and purging the extension tubing last. Unused ports should always be closed, even if there are bidirectional valves. Always check that the stopper caps are properly screwed on. If the devices are connected to a central line, position the manifold at the height of the patient's ear or below, but never above this height. Make a safety loop: part of the extension tubing connecting the device to the catheter must loop below the level of the patient's right ear. This is to avoid the risk of air embolism during or at the end of partial or full gravity infusion. This precautionary measure is mandatory if the manifold cannot be placed below the patient's right ear. Before connection to the patient, check that all ports are closed or connected to another device.

- Connection to the infusion line via luer lock:

- Connection via simple adjustable ring: 1) Move back the sliding ring. 2) Connect the female connector. 3) Move the ring toward the female connector. 4) Secure by screwing.
- Connection via rotating adjustable ring (CAIR SWING): connect by screwing. This connection that secures two infusion devices allows them to rotate freely in relation to each other, thus securing the combination of the two devices.

- Case of devices with Neutraclear™ / Cair Drive bidirectional valve:

- Activation: the activation of the bidirectional valve consists of connecting this valve to a male Luer end fitting, thereby enabling the fluid to flow between the two devices, within a secure and watertight system. Before connecting another male luer device, disinfect the valve septum by rubbing with a pad soaked in alcohol-based disinfectant for 15 seconds. The Neutraclear™ / Cair Drive bidirectional valve is compatible with all Luer slip and Luer lock fittings with a diameter greater than 1.5 mm. To connect a simple Luer, press and turn 1/4 turn to tighten the connection; do not leave unattended. To connect a luer lock, screw all the way on. The maximum number of functional activations for the Neutraclear™ is 600.
- Disconnection: To disconnect, unscrew the male luer end fitting connected to the valve.
- Rinsing: Flush the valve after each injection or sample collection. Flushing is done with a neutral isotonic solution with needle, preferably 10 ml. Flushing must be done in fits and starts, not in a continuous injection. The ideal rinsing volume is between 2.5 and 10 ml. It must be adjusted to the patient. Neutraclear™ valves are transparent and enable the flushing volume to be adjusted through a direct view of its effectiveness.

Neutraclear™ bidirectional valve

Cair Drive bidirectional valve

- Case of devices with C-Flow valve:

- Activation: Before each activation, rub the C-Flow septum with a 70% alcohol-based antiseptic for 15 seconds and leave it to dry.
Connecting a device with a male Luer:
 - If the Luer is a simple Luer: push it in fully then turn 1/4 turn to secure the connection.
 - If the Luer has an adjustable ring, screw it in fully but not too tight.
- Rinsing: Flush the valve after each injection or sample collection referring to hospital protocols for the execution of this step.
- Disconnection: To disconnect, unscrew the male luer end fitting connected to the valve.

C-Flow bidirectional valve

- Duration of use: The protocols in force in hospitals prevail over the durations of use indicated below. Cair guarantees the devices' functional integrity for the following periods:

- 24h if a contrast medium or medicinal products derived from blood are used, or in the presence of a 1.2 µm filter for medialipids
- 96h in the presence of a 0.2µm filter or if presence of a C-Flow valve
- 7 days in all other cases.

- NOTE: Connections must be checked at least once every 24 hours and always before/after every mobilisation of the patient (lifting, positioning in chair, etc.).

5. Caution, precautions for use and contraindication:

Caution: ramp delivered with non-waterproof colored cap, throw it away after connection.

Caution: place the ramp at the maximum height of the patient's right heart atrium.

When connecting two stopcocks or satellite stopcocks to an infusion manifold, extreme care must be taken when handling the stopcocks to prevent the two devices from uncoupling. Likewise, this kind of connection must only be used for patients who are bedbound and cannot be lifted (placed in a chair) and the connection must be checked before and after making their beds.

In all cases, Cair LGL recommends using either a one-piece multi-stopcock manifold, or using stopcocks with a mobile rotating luer connection.

Never connect a male luer cap or use a worn needle or cannula with bidirectional valves as this could damage the device. Should this occur, immediately replace the device.

Glass syringes are not compatible with Neutraclear valves (luer internal diameter < 1.5mm)

Never use pliers. Excessive tightening of connections can damage the luers and make disconnections difficult.

The Neutraclear™, Cair Drive and C-Flow valves are not suitable for use with an undiluted Busulfan solution or a Treanda® solution.

Check any potential incompatibilities and contraindications from drug solution manufacturers before use.

For precaution, close and/or clamp the infusion line when it is not used.

6. Side effects:

The risks incurred by improper use of these devices are risks of embolism or infection.

In the event of a serious incident with the medical device, the user must inform the manufacturer and the relevant health authority of the country where the incident occurred.

It is imperative to keep the incriminated medical devices in order to allow an expertise.

7. Storage/Handling/Disposal:

- Shelf life: Refer to the shelf life described on the label of each device.

- Storage: In a clean, dry, dark area. We recommend to store the device in a place where the temperature is between +5 and +40°C.

- Disposal: according to the institutional procedure for infectious waste and/or waste soiled by biological or chemical liquids.

8. Single use

This device is designed for single use only.

Do not re-sterilize, destroy after use.

CAIR LGL did not design this device to be reprocessed or reused. CAIR LGL therefore cannot certify that reprocessing (e.g.: cleaning or sterilisation) will maintain the device's structural integrity and guarantee the safety of the patient and/or user.

9. Symbols

	Reference number		Batch code		Expiry date
	Sterilised with ethylene oxide		Do not resterilise		Do not reuse. Intended only for single use
	Manufacturer		Always read the instructions for use		Read the instructions for use relating to warning information
	Maximum duration of use (hours)		Maximum duration of use (days)		Pyrogen-free
	Liquid filter with pore diameter		On medical prescription only		Do not use if the packaging is damaged
	Residual volume		Maximum pressure		Manufacturing date
	Temperature limit		Peel-off direction		Quantity
	Store away from sunlight		Fears humidity		Medical device
	Video				

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