FR	Instructions d'utilisation des prolongateurs pour perfusion en PE/PVC et PVC
EN	Instructions for use of PE/PVC and PVC extension tubing for infusion
ES	Instrucciones de uso de los prolongadores para infusión de PE/PVC y PVC
IT	Istruzioni per l'uso delle prolungher per infusione in PE/PVC e PVC
DE	Bedienungsanleitung für Verlängerungen für die Perfusion aus PE/PVC und PVC
DA	Brugsanvisning til infusionsforlængere i PE/PVC og PVC
SV	Instruktioner för användning förlängningsslangar för i PE/PVC och PVC
PT	Instruções de utilização dos prolongadores para perfusão PE/PVC e PVC
NL	Gebruiksaanwijzingen voor verlengslangen voor infusen van PE/PVC en PVC
SL	Navodila za uporabo podaljškov za infundiranje iz PE/PVC in PVC
TR	PE/PVC ve PVC infüzyon amaçlı uzatma hatları kullanım talimatları
FI	Infuusiojatkoliittimien käyttöohjeet, materiaali PE/PVC ja PVC
HR	Upute za uporabu nastavka za infuziju od PE/PVC-a i PVC-a
HU	PE/PVC és PVC perfúzor-hosszabbítók használati utasítása
AR	تعليمات استخدام وصلات إطالة للحقن مصنوعة من البولي إيثيلين/البولي فينيل كلوريد ومن البولي فينيل كلوريد (الرقم
EL	Οδηγίες χρήσης επεκτάσεων για έγχυση από PE/PVC και PVC
NO	Bruerveiledning for forlengelsessett for infusjon i PE/PVC og PVC
PL	Sposób użycia przedłużaczy do zestawów infuzyjnych z PE/PVC i PVC
CS	Pokyny k použití prodlužovacích hadiček pro infúzi z PE/PVC a PVC

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
9. Symboles

1. Indications/Use claims

Extension tubing for one-way or multi-way intravenous infusion is a sterile, single-use medical device. It can include non-return valves, bidirectional valves and/or infusion filters.

Extension tubing is used for the administration of drug substances or sample collections from veins or solely for sample collections from arteries in combination with other infusion devices: infusion pumps, manifolds, stopcocks, bidirectional valves, etc.

These devices are intended to be used by health staff.

These devices are intended to be used on pediatric and adult populations.

2. Composition

Extension tubing can be equipped with a simple luer (the luer end fitting is not screwed onto the device to which it is connected); with a fixed ring (the luer end fitting is screwed onto the device to which it is connected, with no possible movement); or with an adjustable ring (the luer end fitting is screwed onto the device to which it is connected and both devices can rotate in relation to each other).

The devices do not contain DEHP (DOP), latex, BBP, DBP, DIBP or DHP. Any presence of these compounds is a trace amount and is compliant with REACH regulations.

3. Clinical benefit/Performance/Action mechanism

Extension tubing for intravenous infusion helps administer the medicinal product by carrying the drug substances into the patient's bloodstream via a catheter. Extension tubing is also used to collect blood.

The maximum pressure for use of each device is indicated on its technical data sheet.

4. Prerequisites before use and instruction for use

Comply with applicable hygiene rules.

Check the integrity and the cleanliness of the packaging and the device before use; if defective, do not use the device.

Carefully open the blister, taking into account the direction of peeling. Check the expiry date before use; if expired, do not use.

Check all connections (stopper cap) before using the device. Purge the entire device (all ports on the PY) before connecting to the patient.

- Connection to the infusion line via luer lock:

Connection via simple adjustable ring: 1) Move back the sliding ring. 2) Connect the female connector. 3) Move the ring toward the female connector.

4) Secure by screwing.

- Connection via rotating adjustable ring (CAIR SWING): connect by screwing. This connection that secures two infusion devices allows them to rotate freely in relation to each other, thus securing the combination of the two devices.

- Case of extension tubing with non-return valve or anti-siphon valve: non-return valves are open systems, you must therefore clamp downstream of the valve before disconnecting or connecting another device. Risk of entry of air, in particular in central line.

- Case of extension tubing with Neutraclear™ bidirectional valve:

- Activation: the activation of the bidirectional valve consists of connecting this valve to a male luer end fitting, thereby enabling the fluid to flow between the two devices, within a secure and watertight system. Before connecting another male Luer device, disinfect the valve septum by rubbing with a pad soaked in alcohol-based disinfectant for 15 seconds. The Neutraclear™ bidirectional valve is compatible with all Luer slip and Luer lock fittings with an inner diameter greater than 1.5 mm. To connect a simple Luer, press and turn 1/4 turn to tighten the connection; do not leave unattended. To connect a luer lock, screw all the way on. The maximum number of activations for the Neutraclear is 600.
- Disconnection: To disconnect, unscrew the male luer end fitting connected to the valve.
- Rinsing: Flush the valve after each injection or sample collection. Flushing is done with a neutral isotonic solution with needle, preferably 10 ml. Flushing must be done in fits and starts, not in a continuous injection. The ideal rinsing volume is between 2.5 and 10 ml. It must be adjusted to the patient. Neutraclear™ valves are transparent and allow this adjustment of flushing volume through visualisation of the internal cannula.

- Duration of use: The protocols in force in hospitals prevail over the durations of use indicated below. Cair guarantees the devices' functional integrity for the following periods:

- 24h if a contrast medium or medicinal products derived from blood are used, or in the presence of a 1.2 µm filter for medilipids
- 96h in the presence of a 0.2µm filter
- 7 days in all other cases.

- Verification: Connections must be checked at least once every 24 h and before/after each mobilization of the patient (lifting, positioning in chair, etc.).

5. Caution, precautions for use and contraindication

Never connect a male Luer cap or use a worn needle or cannula with bidirectional valves as this could damage the device. Should this occur, immediately replace the device. Glass syringes are not compatible with the valves.

Never use pliers. Excessive tightening of connections can damage the Luers and make connections difficult.

Check any potential incompatibilities and contraindications from drug solution manufacturers before use.

6. Risks and adverse events

The risks incurred by the use of these devices are risks of embolism or infection.

In the event of a serious incident with the medical device, the user must inform the manufacturer and the relevant health authority of the country where the incident occurred.

It is imperative to keep the incriminated medical devices in order to allow an expertise.

7. Storage/Handling/Disposal

Shelf life: refer to the product label. The devices must have a shelf life of 35 months or 59 months.

Storage: in a clean, dry, away from sunlight. We recommend to store the device in a place where the temperature is between +5°C and +40°C.

Waste disposal: To dispose of this device, please refer to the regulation applicable to the disposal of hazardous waste, biological waste or infectious medical waste.


8. Single use

This is a single-use medical device.

Do not resterilize, destroy it after use.

CAIR LGL did not design this device for reprocessing and therefore cannot certify that a reprocessing process maintains the performance and safety characteristics of the device.

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