

Oral care kits



Oral care kit #3

Ref. 401048

- 1 kidney-shaped **basin**
- 4 non-woven **balls** 3 cm, 30 g
- 6 gauze **balls** 2.5x3.5 cm, 17 threads
- 2 **bowls** 60 ml, graduated
- 1 wooden **tongue depressor** adult size
- 1 Kocher plastic **clamp** 14.5 cm, blue
- 4 non-woven **hand towel** 20x20 cm

Oral care accessories

Reference	Description	Quantity per carton
080503	Lemon stick, per sachet of 3	750
080507	Plastic applicator 15x0.4 cm with cotton tip	1000
080509	Star shaped plastic applicator 150x4 mm	1000
500201	Wooden tongue depressor, adult size	5000

Oral care kit #4

Ref. 400076

- 1 **blister pack** 15.5x19.5 cm / 4 cells
- 5 non-woven **swabs** 5x5 cm, 4-ply, 30 g
- 6 non-woven **balls** 3 cm, 30 g
- 1 Kocher plastic **forceps** 14.5 cm, blue
- 1 plastic anatomic **forceps** 12.5 cm green
- 1 adult size wooden **tongue depressor**



product catalogue

Urinary probing Dressings Oral care



Urinary probing kits

Peripheral urinary probing kit

Ref. 405048

- 1 drape 70x90 cm
- 1 sterile urine **bag**, phthalate-free
- 1 **bowl** graduated, 60 ml
- 7 non-woven **balls** 3 cm, 30 g
- 4 non-woven **swabs** 10x10 cm, 4-ply, 30 g
- 2 non-woven **swabs** 10x20 cm, 4-ply, 40 g
- 1 non-woven **swab** 10x20 cm, 4-ply, 40 g
- 1 plastic anatomic **forceps** 12.5 cm blue
- 1 eccentric Luer Lock **syringe**, 10 ml
- 1 NW/PE **drape** 70x70 cm, perforated 6x10 cm, centred
- 1 plastic airlaid **washglove**



Urinary probing kit with washing part #1

Ref. 400108

Washing part

- 1 plastic airlaid **washglove**
- 2 latex exam **gloves** size 7/8, unpowdered, upturned cuffs
- 10 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g

Probing part

- 1 **drape** 50x70 cm
- 1 **blister** pack 4 compartments
- 10 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 1 non-woven **swab** 5x5 cm, 4-ply, 30 g
- 1 **ampoule** of sterile distilled water, 20 ml
- 1 eccentric Luer Lock **syringe**, 20 ml
- 1 pink **needle** 18G 40 mm
- 1 NW/PE **drape** 75x90 cm, with 8 cm perforation at 1/3 + split
- 1 lubricating **gel**
- 1 plastic anatomic **forceps** 12.5 cm
- 1 **bowl** 60 ml, graduated



Urinary probing kit #5

Ref. 400871


- 1 drape 70x90 cm
- 1 NW/PE **drape** 75x90 cm, with 8 cm perforations at 1/3 + split
- 1 **tray** 32.5x17.5x4.5 cm
- 2 **bowls** 150 ml, graduated
- 1 **bowl** 90 ml, graduated
- 15 non-woven **swabs** 10x10 cm, 4-ply, 30 g
- 1 sachet lubricating **gel** 5 g
- 1 eccentric Luer Lock **syringe** 20 ml
- 1 **ampoule** of sterile distilled water 20 ml
- 1 non-woven **hand towel** 30x40 cm



Paper/PE drape,
75x90 cm, perforated
10 cm to 1/3 and pre-cut

Ref. 081307

■ Wadding / PE ■ Opening Pre-cut




Skincare & dressing change kits

With drape

Dressing kit #8

Ref. 400508

- 1 **blister** pack 15.5x19.5 cm / 4 cells
- 5 gauze **balls** 3 cm, 30 g
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 2 plastic anatomic **forceps** 12.5 cm
- 1 paper/PE **drape** 38x45 cm



Dressing kit #25

Ref. 400613

- 1 **blister** pack 15.5x19.5 cm / 4 cells
- 1 paper/PE **drape** 38x45 cm
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 5 gauze **balls** 3 cm, 20 threads
- 1 plastic dressing **clamp** 12.5 cm
- 1 Kocher plastic **clamp** 14.5 cm



Chronic wound deterrentive cleaning kit #1

Ref. 400085

- 1 **blister** pack 15.5x19.5 cm / 4 cells
- 1 paper/PE **drape** 38 x 45 cm
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 1 Kocher plastic **clamp** 14.5 cm, blue
- 1 plastic dressing **clamp** 12.5 cm
- 1 metal **scissors** plastic loops, 12 cm
- 1 plastic **stylet** 14 cm
- 1 plastic **curette** 13.5 cm



Without drape

Dressing kit #2

Ref. 400502

- 1 **blister** pack 15.5x19.5 cm / 4 cells
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 5 gauze **balls** 3 cm, 20 threads
- 2 plastic anatomic **forceps** 12.5 cm

