

Infusion kit for implantable port



Ref. 400029-1

- 1 tray
- 1 non-woven **gown** with jersey cuffs
- 4 nitrile **gloves** size 7/8, unpowdered, upturned cuffs
- 2 type IIR elastic surgical **masks**
- 1 **drape** 50 x 75 cm 11 cm perforated to 1/3 adhesive + pre-cut
- 15 non-woven **swabs** 7.5 x 7.5 cm, 4-ply, 30 g
- 1 Luer Lock centred **syringe** 20 ml
- 1 Luer Lock 3-way extension **valve** 30 cm
- 1 Luer Lock 3-way male female **valve**
- 1 **ampoule** NaCl 0.9% 20 ml
- 2 double **strips** 15 x 1.25 cm
- 1 safety transfer **needle** 18G 40 mm
- 1 transparent **adhesive film dressing** 12 x 14 cm

Single sterile range

Reference	Description	Sale unit	Quantity per carton
062618	Double strips	Box of 20 sachets of 2	1000
370107	Foam sponge stick, 15 cm	Sachet of 1	90
401107	Drape, 37 x 45 cm	Per box	270
350112	Scissors, 12 cm	Bag of 1	400
800100	Latex-free tourniquet	Per unit	2500
800116	Sterile latex-free canvas flat tourniquet	Bag of 1	110
360101	Plastic anatomic forceps, 12.5 cm, blue	Bag of 1	120
360132	Plastic clamping forceps, 12 cm	Bag of 1	240
350202	Tray, 26.5 x 16.2 x 4.5 cm	Bag of 1	40
350203	Tray, 32.5 x 17.5 x 4.5 cm	Bag of 1	34
390139	Pre-filled syringe, NaCl 0.9% 5 ml, sterile (class IIa)	Per unit	360
390129	Pre-filled syringe, NaCl 0.9%, 10 ml (class IIa)	Per unit	320
390135	Pre-filled syringe, NaCl 0.9%, 20 ml (class IIa)	Per unit	150
390136	Pre-filled syringe, NaCl 0.9%, 10 ml (class III)	Per unit	240
390325	Sterile neutral pressure bidirectional valve	Per unit	200
060921	Statlock catheter maintenance system	Bag of 1	50



product catalogue

Infusion

Implantable port - Central venous line
- Peripheral venous line

CVL infusion kit #1

Ref. 400049-2

- 1 table **drape** 50x70 cm
- 4 nitrile examination **gloves**, size 7/8, upturned cuff
- 2 type IIR elastic surgical **masks**
- 1 white **charlotte**
- 1 **drape** 70x90 cm, with side adhesive tape 70 cm
- 15 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 2 Luer Lock **syringes**, 20 ml
- 2 safety transfer **needles** 18G 40 mm
- 1 x 15 L biomedical **waste bag**, yellow



Chemotherapy kit #5

Ref. 400414

- 1 blister pack
- 1 **drape** 50x50 cm, folded
- 1 double **strip** 15x1.25 cm
- 1 adhesive dressing with central **swab** 10x15 cm
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 5 non-woven **balls** 3 cm, 30 g
- 1 Kocher plastic **clamp** 14.5 cm
- 1 plastic anatomic **forceps** 12.5 cm
- 1 **drape** 50x50 cm, perforated 11 cm eccentric + 3 adhesives + pre-cut



Infusion kit for implantable port #8

Ref. 400266-1

- 1 table **drape** 50x70 cm
- 1 green non-woven **gown** with jersey cuffs, 30 g
- 4 nitrile examination **gloves**, size 7/8, upturned cuff
- 2 type IIR surgical **masks**
- 1 **drape** 50x50 cm perforated 11 cm eccentric + 3 adhesives + pre-cut
- 15 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 1 Luer Lock 3-way **extension valve**, 30 cm, mobile
- 1 Luer Lock **male female** 3-way **valve**
- 1 double **strip** 15x1.25 cm
- 1 transparent **adhesive film** dressing 12x14 cm
- 1 pre-filled NaCl **syringe**, 10 ml



Dressing change for PICC line #2

Ref. 400679-1

- 1 pre-filled **syringe**, NaCl 0.9%, 20 ml
- 2 type IIR elastic surgical **masks**
- 1 table **drape** 50x50 cm
- 5 non-woven **swabs** 7.5x7.5 cm, 4-ply, 30 g
- 12 non-woven **balls** 3 cm, 30 g
- 1 plastic anatomic **forceps** 12.5 cm, green
- 1 plastic anatomic **forceps** 12.5 cm, blue
- 1 double **strip** 15x1.25 cm
- 1 transparent **adhesive film** dressing 10x12 cm

