

Clamping forceps

Unscrew cap


An innovative system

Thanks to **the notches**, the clamping forceps can also be used **to unscrew plugs or valves** during handling.

Soon available in all our sets containing clamping forceps..



Unscrew cap



REF 360112

Packaging :
available in sets.

Scan & discover
the video !



07-2023

This is a **Class I medical device**.

This medical device is a **regulated health product** which carries, under this regulation, **the CE marking** which RAFFIN MEDICAL is the manufacturer.



This clamp is **used to close and/or regulate the flow of fluids** in dialysis and intravenous infusion.
Can be used **for all types of flexible plastic tubing** that need to be closed.

