

NEW RANGE

Elastic tubular bandages

Application areas

Finger, hand, wrist, foot, arm, head, abdomen, chest...

Available sizes


0.5 to 11 (adult/child)

Box content

1 roll of 6 m (unstretched)
(25 m stretched)

Composition

Polyester: 40 %
Polyamide: 20 %
Spandex: 40 %



The elastic tubular bandage is **designed to fit all areas of the body**. Lightweight, flexible, and stretchable, it ensures effective retention of dressings.

Thanks to its innovative structure, it provides optimal support and high comfort for the patient. It can be cut at every length to meet specific needs.



Reference	Size	Unstretched width	Stretched width	Intended use
020623	0,5	1 cm	7 cm	Finger
020624	1	1,4 cm	15,5 cm	Child's hand
020625	2	1,6 cm	24 cm	Child's wrist/foot
020626	3	2 cm	27,5 cm	Hand
020627	4	2,2 cm	31 cm	Hand/forearm/arm
020628	5	2,7 cm	37 cm	Child's leg/head
020629	6	3,9 cm	51 cm	Head Abdomen/chest S
020630	7	5,2 cm	70 cm	Abdomen/chest M
020631	8	6,2 cm	81 cm	Abdomen/chest L
020632	9	7,2 cm	90 cm	Abdomen/chest XL
020633	10	8,1 cm	105 cm	Abdomen/chest XXL
020634	11	10,2 cm	126 cm	Abdomen/chest XXXL

This Class I medical device (non-sterile) is a regulated health product bearing the **CE** marking, for which RAFFIN MEDICAL is the European authorized representative and importer.

